Apparatus: Electrical muffle furnace “economic chamber furnace”

Models: HD 230

Execution: PAD

Feeding tension: 220v monophasic

Consumption: According to characteristic plate placed on each furnace model

GENERAL RULES:

- A)** To work within the temperature regimes assigned to the furnace and not to exceed the established limit temperatures, determined by the disposition of the heating elements and the required safety margins.

Furnace model **HD 230 “PAD”**

1.200°C Maximum temperature.

- B)** The fact that an electric furnace is offered up to determined maximum work temperatures does not mean that it cannot exceed them if it is not under supervision or automatic control, which can lead to the fusion of the heating elements. The limit work temperatures assigned correspond to the maximum admissible based on the heating elements within the normal functioning guarantees and the suitable safety margins. Considering that there is always a temperature differential between the registered temperature in the useful chamber and that of the heating element, always being superior the latter, in order to compensate the furnace heat losses and to increase relatively quickly the chamber temperature that require the work needs, we advise that, if the maximum temperatures that the furnace can reach are not necessary, it is better not to use them without a determined purpose as this always causes a reduction of the normal life of the heating element due to its work at the most critical points.

In the furnaces with resistances that work with direct radiation, this basic characteristic should not be modified. The obstruction of the normal radiation of the heating element with compact objects not specially designed that produce an excessive heat concentration in the heating elements, means overtemperature and fusion risks. In general it should be avoided all those thermal processes that chemically or physically attack the heating resistances and in any case it will be avoided the direct contact with strange materials to the ceramic parts of the muffle or supports.

- C)** To dispose of the suitable temperature control tools to know at any time the exact chamber temperature as well as to apply the suitable regulation apparatus when it is necessary to carry out thermal processes that require to stabilise the furnace at high temperatures for a long time. If due to long furnace use or other causes there are doubts on the marking accuracy of the pyrometric equipment, it is advisable to contrast it with a known fidelity pattern.
- D)** We always recommend the use of original spare parts in the HOBERSAL furnaces and the reparation services of the manufacturer himself; nobody better than him can leave the furnace again with the same conditions and guarantees of first use.

The fulfilment of such rules and the careful treatment of the apparatus will allow to get from the furnace a suitable performance, long life and elimination of annoying breakdowns.

As manufacturers, we recommend the users that in the interest of both, they grant these rules the necessary attention and importance. Thanks.

INSTRUCTIONS

Connection: It will be connected to the supply system at 220v between active phases. We recommend to establish the ground connection that the pin incorporates – origin plug. This ground connection has to be always a real ground, never the neutral of the mains. In the event that there were no ground installation, it is preferable not to connect it as it is not essential for the normal functioning of the apparatus, although it is very recommendable as at high temperatures heat-resistant materials loose partially their dielectric value and in some cases there can be a little electrical derivation to the general chassis of the apparatus, a leak originated through the heat-resistant parts that support the electrical resistances and that are subject to maximum temperature.

This normal phenomenon is already planned in the furnace conception, with the due safety isolation, in the electrical aspect in the parts of the apparatus manipulated by the operator, and it is totally eliminated establishing the regulated and correct ground connection.

Chimney: At the back of the furnace, there is a round assembly destined to the fixation of the self-sucking chimney supplied as an accessory. If this accessory has not been ordered and the furnace does not have any chimney, the hole can be left uncovered. If you wish, you can cover it with ceramic fibre or another fireproof material, but we repeat that it is not necessary to do it for the normal functioning of the furnace.

BOTTOM TRAY: It protects the plate and/or resistance of the furnace base against product spills or electrical contacts. It should always be placed with the ribs in the lower part (on the furnace bottom) to allow a freer heat radiation of the heating plate of the furnace base.

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