

GENERAL INFORMATION

Product name : Beaker, low form, PMP

Description : Made of Polymethylpentene (PMP/TPX®). Higher transparency, chemical resistance and better heat tolerance than PP. Autoclavable. Easy to read silk screened graduations. Stackable. According to ISO 7056:1981

TECHNICAL DATA

reference	capacity	graduation	Øouter (mm)	h (mm)	pcs/pack
BKLX-050-012	50 ml	5 ml	49.4	58.4	12
BKLX-100-012	100 ml	10 ml	63.2	72.5	12
BKLX-250-006	250 ml	25 ml	83.1	92.1	6
BKLX-600-004	600 ml	50 ml	110.3	131	4
BKLX-1K0-003	1000 ml	50 ml	129.4	149.2	3
BKLX-2K0-001	2000 ml	100 ml	158.5	190.2	1
BKLX-5K0-001	5000 ml	500 ml	202.9	260.5	1
BKLX-10K-001	10000 ml	1000 ml	255.4	326.6	1

PACKAGING AND LOGISTICS

reference	vol (l)	kg	TARIC	GTIN
BKLX-050-012	0,01	0,01	39233010	08434868007499
BKLX-100-012	0,01	0,01	39233010	08434868007505
BKLX-250-006	0,01	0,01	39233010	08434868007529
BKLX-600-004	0,01	0,01	39233010	08434868007550
BKLX-1K0-003	0,01	0,01	39233010	08434868007512
BKLX-2K0-001	0,01	0,01	39233010	08434868007536
BKLX-5K0-001	0,01	0,01	39233010	08434868007543
BKLX-10K-001	12	1	39233010	08434868086302

PRODUCT PHOTO



INFORMACIÓN GENERAL

Nombre del producto : Vaso de precipitados forma baja de PMP**Descripción :** Fabricado en Polimetilpenteno (PMP/TPX®). Mayor transparencia, resistencia química y tolerancia al calor que el PP. Autoclavable. Graduación serigrafiada y fácil de leer. Apilable. Conforme a ISO 7056:1981

DATOS TÉCNICOS

referencia	capacidad	graduación	Øexterior(mm)	h (mm)	unidades por ref.
BKLX-050-012	50 ml	5 ml	49.4	58.4	12
BKLX-100-012	100 ml	10 ml	63.2	72.5	12
BKLX-250-006	250 ml	25 ml	83.1	92.1	6
BKLX-600-004	600 ml	50 ml	110.3	131	4
BKLX-1K0-003	1000 ml	50 ml	129.4	149.2	3
BKLX-2K0-001	2000 ml	100 ml	158.5	190.2	1
BKLX-5K0-001	5000 ml	500 ml	202.9	260.5	1
BKLX-10K-001	10000 ml	1000 ml	255.4	326.6	1

EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
BKLX-050-012	0,01	0,01	39233010	08434868007499
BKLX-100-012	0,01	0,01	39233010	08434868007505
BKLX-250-006	0,01	0,01	39233010	08434868007529
BKLX-600-004	0,01	0,01	39233010	08434868007550
BKLX-1K0-003	0,01	0,01	39233010	08434868007512
BKLX-2K0-001	0,01	0,01	39233010	08434868007536
BKLX-5K0-001	0,01	0,01	39233010	08434868007543
BKLX-10K-001	12	1	39233010	08434868086302

FOTO DEL PRODUCTO



INFORMATIONS GÉNÉRALES

Nom produit : Bécher forme basse en PMP

Description : Fabriqué en Polyméthylpentène (PMP/TPX®). Transparence, résistance chimique et tolérance à la chaleur plus grandes que le PP. Autoclavable. Graduations en sérigraphie facile à lire. Empilable. Conforme ISO 7056:1981

DONNÉES TECHNIQUES

référence	capacité	graduation	Øextérieur (mm)	h (mm)	unités par ref.
BKLX-050-012	50 ml	5 ml	49.4	58.4	12
BKLX-100-012	100 ml	10 ml	63.2	72.5	12
BKLX-250-006	250 ml	25 ml	83.1	92.1	6
BKLX-600-004	600 ml	50 ml	110.3	131	4
BKLX-1K0-003	1000 ml	50 ml	129.4	149.2	3
BKLX-2K0-001	2000 ml	100 ml	158.5	190.2	1
BKLX-5K0-001	5000 ml	500 ml	202.9	260.5	1
BKLX-10K-001	10000 ml	1000 ml	255.4	326.6	1

EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
BKLX-050-012	0,01	0,01	39233010	08434868007499
BKLX-100-012	0,01	0,01	39233010	08434868007505
BKLX-250-006	0,01	0,01	39233010	08434868007529
BKLX-600-004	0,01	0,01	39233010	08434868007550
BKLX-1K0-003	0,01	0,01	39233010	08434868007512
BKLX-2K0-001	0,01	0,01	39233010	08434868007536
BKLX-5K0-001	0,01	0,01	39233010	08434868007543
BKLX-10K-001	12	1	39233010	08434868086302



INFORMAZIONE GENERALE

Nome del prodotto : Becher formato basso in PMP**Descrizione :** Realizzato in polimetilpentene (PMP/TPX®). Maggiore trasparenza e resistenza chimica e termica al calore rispetto al PP. Autoclavabile. Graduazione serigrafata facilmente leggibile. Impilabile. Conforme a ISO 7056:1981

DATI TECNICI

referenza	capacità	graduazione	Øesterno (mm)	h (mm)	unità per ref.
BKLX-050-012	50 ml	5 ml	49.4	58.4	12
BKLX-100-012	100 ml	10 ml	63.2	72.5	12
BKLX-250-006	250 ml	25 ml	83.1	92.1	6
BKLX-600-004	600 ml	50 ml	110.3	131	4
BKLX-1K0-003	1000 ml	50 ml	129.4	149.2	3
BKLX-2K0-001	2000 ml	100 ml	158.5	190.2	1
BKLX-5K0-001	5000 ml	500 ml	202.9	260.5	1
BKLX-10K-001	10000 ml	1000 ml	255.4	326.6	1

IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
BKLX-050-012	0,01	0,01	39233010	08434868007499
BKLX-100-012	0,01	0,01	39233010	08434868007505
BKLX-250-006	0,01	0,01	39233010	08434868007529
BKLX-600-004	0,01	0,01	39233010	08434868007550
BKLX-1K0-003	0,01	0,01	39233010	08434868007512
BKLX-2K0-001	0,01	0,01	39233010	08434868007536
BKLX-5K0-001	0,01	0,01	39233010	08434868007543
BKLX-10K-001	12	1	39233010	08434868086302

FOTO DEL PRODOTTO

